

GFRA-1 protein (STJP000149)

STJP000149

GENERAL INFORMATION

Product Type Proteins
Short Description GFRA-1-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 0.9 mg/ml
Formulation PBS pH 7.4
Purification Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID 2674
Gene Symbol GFRA1
Uniprot ID GFRA1_HUMAN
Immunogen Sequence