

## GFRA-4 protein (STJP000151)

STJP000151

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** GFRA-4-protein was developed for use in research applications.  
**Host/Source**

### PRODUCT PROPERTIES

**Concentration** 1.3 mg/ml  
**Formulation** PBS pH 7.4  
**Purification** Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media  
**Dilution Range**  
**Storage Instruction** Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

**Gene ID** 64096  
**Gene Symbol** GFRA4  
**Uniprot ID** GFRA4\_HUMAN  
**Immunogen Sequence**