

RET-His protein (STJP000152)

STJP000152

GENERAL INFORMATION

Product Type Proteins
Short Description RET-His-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 0.23 mg/ml
Formulation PBS pH 7.4
Purification Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID 5979
Gene Symbol RET
Uniprot ID RET_HUMAN
Immunogen Sequence