

Cardiac Troponin complex protein (STJP000154)

STJP000154

GENERAL INFORMATION

Product Type Proteins
Short Description Cardiac Troponin complex-protein was developed from cho cells. For use in research applications.
Host/Source CHO cells

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation 25 mM Tris, pH 7.5, 500 mM NaCl 5 mM CaCl₂, 15 mM Beta-mercaptoethanol
Purification Antibody-affinity purification following desalting
Dilution Range >90%
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [7137](#)
Gene Symbol [TNNI3](#)
Uniprot ID [TNNI3_HUMAN](#)
Immunogen Sequence