

ACE2 protein (hFc) (STJP000155)

STJP000155

GENERAL INFORMATION

Product Type	Proteins
Short Description	ACE2-hFc protein was developed from human and has a target region of hFc. For use in research applications.
Host/Source	Human

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	1. PBS pH 7.4, 2. 20 mM Na-citrate, 150 mM NaCl; pH 6.0
Purification	Purified by Protein A affinity chromatography following gel filtration. Protein is sterile-filtered through 0.22 µm filter.
Dilution Range	> 95%
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	59272
Gene Symbol	ACE2
Uniprot ID	ACE2_HUMAN
Immunogen Sequence	