

## ACE2 protein (His-Tag) (STJP000156)

STJP000156

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** ACE2-His-Tag protein was developed from human and has a target region of His-Tag. For use in research applications.  
**Host/Source** Human

### PRODUCT PROPERTIES

**Concentration** 1 mg/ml  
**Formulation** PBS pH 7.4  
**Purification** Metal-affinity chromatography following gel filtration. Protein is sterile-filtered through 0.22 µm filter.  
**Dilution Range** > 80%  
**Storage Instruction** Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

**Gene ID** [59272](#)  
**Gene Symbol** [ACE2](#)  
**Uniprot ID** [ACE2\\_HUMAN](#)  
**Immunogen Sequence**