

Recombinant SARS-CoV-2 Nucleocapsid protein (His-Tag) (STJP000158)

STJP000158

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant SARS-CoV-2 Nucleocapsid-His-Tag protein was developed from sars-cov-2 and has a target region of His-Tag. For use in research applications.
Host/Source	SARS-CoV-2

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Metal-affinity chromatography following gel filtration. Protein is sterile-filtered through 0.22 µm filter.
Dilution Range	>90%
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence