

## Recombinant murine leukemia virus p30 protein (STJP000160)

STJP000160

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant murine leukemia virus p30-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	0.25 mg/ml
<b>Formulation</b>	PBS, 300 mM NaCl pH 7.4
<b>Purification</b>	Purified by Ni affinity chromatography following gel-filtration
<b>Dilution Range</b>	
<b>Storage Instruction</b>	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen Sequence</b>	