

Artemin protein (STJP000161)

STJP000161

GENERAL INFORMATION

Product Type Proteins
Short Description Artemin-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation 10 mM Na-citrate; 150 mM NaCl pH 4.5
Purification Purified from serum-free CHO growth media by ion-exchange chromatography and gel-filtration. Protein is sterile-filtered through 0.22 µm filter
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [9048](#)
Gene Symbol [ARTN](#)
Uniprot ID [ARTN_HUMAN](#)
Immunogen
Sequence