

RNase 7 protein (STJP000162)

STJP000162

GENERAL INFORMATION

Product Type Proteins
Short Description RNase 7-protein was developed from cho cells. For use in research applications.
Host/Source CHO cells

PRODUCT PROPERTIES

Concentration 1.9 mg/ml
Formulation 10 mM Na-acetate, pH 5.0; 150 mM NaCl
Purification Purified by affinity chromatography and gel-filtration from serum-free CHO growth media.
Dilution Range
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID [84659](#)
Gene Symbol [RNASE7](#)
Uniprot ID [RNASE7_HUMAN](#)
Immunogen Sequence