

Recombinant SARS-CoV-2 S1 RBD2 protein (STJP000163)

STJP000163

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant SARS-CoV-2 S1 RBD2-protein was developed for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation PBS pH 7.4
Purification Metal-affinity chromatography following gel filtration. Protein is sterile-filtered through 0.22 µm filter.
Dilution Range > 95%
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence