

Recombinant SARS-CoV-2 Spike S2 ECD-His protein (STJP000165)

STJP000165

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant SARS-CoV-2 Spike S2 ECD-His-protein was developed from sars-cov-2. For use in research applications.
Host/Source	SARS-CoV-2

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter.
Dilution Range	> 90%
Storage	Store at -70°C upon receipt. It is strongly recommended to avoid repeated freeze-thaw cycles. Aliquot the protein, freeze quickly in liquid nitrogen and thaw only once for your experiments.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence