

## Recombinant SARS-CoV-2 Trimeric Spike protein (STJP000166)

STJP000166

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** Recombinant SARS-CoV-2 Trimeric Spike-protein was developed for use in research applications.  
**Host/Source**

### PRODUCT PROPERTIES

**Concentration** 1 mg/ml  
**Formulation** PBS pH 7.4  
**Purification** Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated  
**Dilution Range** >95%  
**Storage Instruction** Store at -70°C. Avoid multiple freeze-thaw cycles

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen Sequence**

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081