

SARS-CoV-2 Spike RBD2 Indian B.1.617.1 protein (Kappa) (STJP000170)

STJP000170

GENERAL INFORMATION

Product Type	Proteins
Short Description	SARS-CoV-2 Spike RBD2 Indian B.1.617.1-Kappa protein was developed from cho cells and has a target region of Kappa. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter.
Dilution Range	> 95%
Storage	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence