

SARS-CoV-2 Spike S1 Indian B.1.617.2 protein (Delta) (STJP000171)

STJP000171

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | SARS-CoV-2 Spike S1 Indian B.1.617.2-Delta protein was developed from cho cells and has a target region of Delta. For use in research applications. |
| Host/Source | CHO cells |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | 1 mg/ml |
| Formulation | PBS pH 7.4 |
| Purification | Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter. |
| Dilution Range | >95% |
| Storage Instruction | Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence