

## SARS-CoV-2 Spike S1 Indian B.1.617.1 protein (Kappa) (STJP000172)

STJP000172

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	SARS-CoV-2 Spike S1 Indian B.1.617.1-Kappa protein was developed from cho cells and has a target region of Kappa. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	1 mg/ml
<b>Formulation</b>	PBS pH 7.4
<b>Purification</b>	Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter.
<b>Dilution Range</b>	>95%
<b>Storage</b>	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Sequence</b>