

SARS-CoV-2 Spike S1 Brazil P1 protein (Gamma) (STJP000174)

STJP000174

GENERAL INFORMATION

Product Type	Proteins
Short Description	SARS-CoV-2 Spike S1 Brazil P1-Gamma protein was developed from cho cells and has a target region of Gamma. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter.
Dilution Range	>95%
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence