

Recombinant Interleukin 6 protein protein (STJP000183)

STJP000183

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant Interleukin 6 protein-protein was developed from cho cells. For use in research applications.
Host/Source CHO cells

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation PBS pH 7.4
Purification
Dilution Range > 95%
Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence