

SARS-CoV-2 Trimeric Spike Californian B.1.429.1 protein (STJP000185)

STJP000185

GENERAL INFORMATION

Product Type Proteins
Short Description SARS-CoV-2 Trimeric Spike Californian B.1.429.1-protein was developed from cho cells. For use in research applications.
Host/Source CHO cells

PRODUCT PROPERTIES

Concentration 1 mg/ml
Formulation PBS pH 7.4
Purification Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated
Dilution Range >95%
Storage Instruction Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence