

Human Angiopoietin-2 protein (Recombinant) {Biotin} (STJP000188)

STJP000188

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Angiopoietin-2-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Angiopoietin-2 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Angiopoietin-2 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	285
Gene Symbol	ANGPT2
Uniprot ID	ANGP2_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081