

Human Interleukin-4 protein (Recombinant) {Biotin} (STJP000191)

STJP000191

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Interleukin-4-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant IL4 was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	3565
Gene Symbol	IL4
Uniprot ID	IL4_HUMAN
Immunogen Sequence	