

Human CD32a (FcγRIIa) protein (Recombinant) {Biotin} (STJP000197)

STJP000197

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD32a (FcγRIIa)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD32a is supplied as a lyophilized 0.2 μg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD32a, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	2212
Gene Symbol	FCGR2A
Uniprot ID	FCG2A_HUMAN
Immunogen Sequence	