

Human RANTES-CCL5 protein (Recombinant) {Biotin} (STJP000200)

STJP000200

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human RANTES-CCL5-protein was developed from e.coli. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Rantes/CCL5was lyophilized from a.2 Mu m filtered PBS solution.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at-20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	6352
Gene Symbol	CCL5
Uniprot ID	CCL5_HUMAN
Immunogen Sequence	