

Human SDF-1 alpha protein (Recombinant) {Biotin} (STJP000208)

STJP000208

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human SDF-1 alpha-protein was developed from e.coli. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant SDF-1 alpha was lyophilized from a.2 Mu m filtered PBS, pH.4.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	6387
Gene Symbol	CXCL12
Uniprot ID	SDF1_HUMAN
Immunogen	
Sequence	