

Human MCP-3 (CCL7) protein (Recombinant) {Biotin} (STJP000209)

STJP000209

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human MCP-3 (CCL7)-protein was developed from e.coli. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	MCP-3 was lyophilized from a 2 µg/ml filtered mM PB solution pH.5, mM NaCl.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	6355
Gene Symbol	CCL8
Uniprot ID	CCL8_HUMAN
Immunogen Sequence	