

Human HER4 (ERBB4) protein (Recombinant) {Biotin} (STJP000211)

STJP000211

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human HER4 (ERBB4)-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human ERBB4 was lyophilized from a 0.2 Mu m filtered PBS pH 7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable los
Instruction	

TARGET INFORMATION

Gene ID	2066
Gene Symbol	ERBB4
Uniprot ID	ERBB4_HUMAN
Immunogen Sequence	