

Human MEC (CCL28) protein (Recombinant) {Biotin} (STJP000217)

STJP000217

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human MEC (CCL28)-protein was developed from e.coli. For use in research applications.
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CCL-28 was lyophilized from PBS pH.4.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56477
Gene Symbol	CCL28
Uniprot ID	CCL28_HUMAN
Immunogen	
Sequence	