

Human Interferon-gamma protein (Recombinant) (STJP000222)

STJP000222

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Interferon-gamma-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Interferon gamma was lyophilized from a 2 µg/ml filtered concentrated (1mg/ml) solution in PBS, pH.7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	3458
Gene Symbol	IFNG
Uniprot ID	IFNG_HUMAN
Immunogen	
Sequence	