

Human Interleukin-1 alpha protein (Recombinant) (STJP000223)

STJP000223

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Interleukin-1 alpha-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Interleukin-1 alpha was lyophilized from a 2 µg/ml filtered PBS solution pH 7.0.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	3552
Gene Symbol	IL1A
Uniprot ID	IL1A_HUMAN
Immunogen	
Sequence	