

Mouse GM-CSF protein (Recombinant) (STJP000224)

STJP000224

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse GM-CSF-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Mouse GM-CSF was lyophilized from a 0.2 µg/ml filtered solution in PBS, pH 7.0.
Purification	
Dilution Range	>96%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss.
Instruction	

TARGET INFORMATION

Gene ID	1437
Gene Symbol	CSF2
Uniprot ID	CSF2_HUMAN
Immunogen Sequence	