

Human GM-CSF protein (Recombinant) (STJP000225)

STJP000225

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant-Human GM-CSF-protein was developed from cho cells. For use in research applications.
Host/Source CHO cells

PRODUCT PROPERTIES

Concentration
Formulation Human GM-CSF was lyophilized from a 2 µg/ml filtered solution in PBS, pH7.0.
Purification
Dilution Range >96%, as determined by SDS-PAGE and HPLC
Storage The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction

TARGET INFORMATION

Gene ID 1437
Gene Symbol CSF2
Uniprot ID CSF2_HUMAN
Immunogen
Sequence