

Mouse LIF protein (Recombinant) (STJP000229)

STJP000229

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse LIF-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Leukemia Inhibitory Factor was lyophilized from a.2 Mu m filtered PBS pH.5.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC.
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	16878
Gene Symbol	Lif
Uniprot ID	LIF_MOUSE
Immunogen	
Sequence	