

Human Oncostatin M protein (Recombinant) (STJP000232)

STJP000232

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Oncostatin M-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Oncostatin M was lyophilized from a.2 Mu m filtered concentrated (1mg/ml) PBS solution, pH7.0.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	5008
Gene Symbol	OSM
Uniprot ID	ONCM_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081