

## Human VEGF 165AA protein (Recombinant) (STJP000234)

STJP000234

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human VEGF 165AA-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant VEGF was lyophilized from 2 µg/ml filtered mM PB solution, pH 7.0.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	7422
<b>Gene Symbol</b>	VEGFA
<b>Uniprot ID</b>	VEGFA_HUMAN
<b>Immunogen</b>	
<b>Sequence</b>	