

Human Stem Cell Factor (mammalian) protein (Recombinant) (STJP000236)

STJP000236

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant-Human Stem Cell Factor (mammalian)-protein was developed from cho cells. For use in research applications.
Host/Source CHO cells

PRODUCT PROPERTIES

Concentration
Formulation Recombinant SCF was lyophilized from a.2 Mu m filtered solution inPBSPH7.0.
Purification
Dilution Range >97%, as determined by SDS-PAGE and HPLC
Storage The lyophilized protein is stable for at least years from date of receipt at-20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months, with a carrier protein without detectable loss
Instruction

TARGET INFORMATION

Gene ID 4254
Gene Symbol KITLG
Uniprot ID SCF_HUMAN
Immunogen
Sequence