

Human PLGF protein (Recombinant) (STJP000238)

STJP000238

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human PLGF-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant PLGF was lyophilized from a 0.2 Mu m filtered solution inPBSpH 7.0.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at-20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months, with a carrier protein without detectable los
Instruction	

TARGET INFORMATION

Gene ID	4254
Gene Symbol	KITLG
Uniprot ID	SCF_HUMAN
Immunogen	
Sequence	