

Mouse VEGF 164AA protein (Recombinant) (STJP000241)

STJP000241

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant-Mouse VEGF 164AA-protein is for use in research applications.
Host/Source

PRODUCT PROPERTIES

Concentration
Formulation
Purification
Dilution Range
Storage Instruction

TARGET INFORMATION

Gene ID [22339](#)
Gene Symbol [Vegfa](#)
Uniprot ID [VEGFA_MOUSE](#)
Immunogen Sequence