

Cynomolgus ACE-2 protein (Recombinant) (Fc-Tag) (STJP000243)

STJP000243

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Cynomolgus ACE-2-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant ACE2 was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	