

Rat IL-20RB protein (Recombinant) (STJP000245)

STJP000245

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Rat IL-20RB-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Rat Interleukin-20 receptor beta is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant rat IL-20RB, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	53832
Gene Symbol	IL20RA
Uniprot ID	I20RA_HUMAN
Immunogen	
Sequence	