

Rat IL-9 receptor (CD129) protein (Recombinant) (STJP000247)

STJP000247

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Rat IL-9 receptor (CD129)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant rat IL-9R (CD129) is lyophilized from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Rat IL9R can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence