

Hepatitis C Virus Nonstructural protein 4 protein (Recombinant) (STJP000250)

STJP000250

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Hepatitis C Virus Nonstructural protein 4-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant BCMA was lyophilized from a 2 µm filtered mm Tris solution pH 5.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this product can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss o
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	