

Human CD3 Epsilon Delta heterodimer protein (Recombinant) (STJP000253)

STJP000253

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD3 Epsilon Delta heterodimer-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human CD3 Epsilon-CD3 Delta is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD3 heterodimer can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	916
Gene Symbol	CD3E
Uniprot ID	CD3E_HUMAN
Immunogen Sequence	