

COV2 Spike protein RBD B.1.1.529 (Omicron) protein (Recombinant) (His-Tag) (STJP000263)

STJP000263

GENERAL INFORMATION




Product Type	Proteins
Short Description	Recombinant-COV2 Spike protein RBD B.1.1.529 (Omicron)-His-Tag protein was developed from hek293 and has a target region of His-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant COVID-19 Spike RBD was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081