

2019 nCoV/COVID19 Spike protein S1 protein (Recombinant) (Fc-Tag) (STJP000268)

STJP000268

GENERAL INFORMATION


Product Type	Proteins
Short Description	Recombinant-2019 nCoV/COVID19 Spike protein S1-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant COVID-19 Spike S1 was lyophilized from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence



1µg loaded Tris-Bis PAGE under reducing condition. The purity is greater than 95%.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081