

## 2019 nCoV/COVID19 Spike protein S1 protein (Recombinant) (His-Tag) (STJP000269)

STJP000269

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-2019 nCoV/COVID19 Spike protein S1-His-Tag protein was developed from hek293 and has a target region of His-Tag. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant COVID-19 Spike S1 was lyophilized from 0.2 µm filtered PBS solution pH7.4
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.
<b>Instruction</b>	

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Sequence**



1µg loaded Tris-Bis PAGE under reducing condition. The purity is greater than 95%.

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081