

Mouse IL-13 alpha 1 receptor (IL13RA1) protein (Recombinant) (STJP000279)

STJP000279

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-13 alpha 1 receptor (IL13RA1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse IL-13 receptor alpha 1 is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant mouse IL13RA1, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16164
Gene Symbol	Il13ra1
Uniprot ID	I13R1_MOUSE
Immunogen	
Sequence	