

Human Tissue factor pathway inhibitor 2 protein (Recombinant) (STJP000287)

STJP000287

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Tissue factor pathway inhibitor 2-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Tissue factor pathway inhibitor 2 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Tissue factor pathway inhibitor can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	21789
Gene Symbol	Tfpi2
Uniprot ID	TFPI2_MOUSE
Immunogen	
Sequence	