

Human KLF4 protein (Recombinant) (STJP000289)

STJP000289

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human KLF4-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human KLF4 was lyophilized from a 0.2 µm filtered PBS solution pH 7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss.
Instruction	

TARGET INFORMATION

Gene ID	9314
Gene Symbol	KLF4
Uniprot ID	KLF4_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081