

Mouse Tissue factor pathway inhibitor protein (Recombinant) (STJP000290)

STJP000290

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse Tissue factor pathway inhibitor-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse Tissue factor pathway inhibitor is lyophilized from 0.2 M unfiltered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Tissue factor pathway inhibitor can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	21788
Gene Symbol	Tfpi
Uniprot ID	TFPI1_MOUSE
Immunogen	
Sequence	