

Mouse IL-31R beta (OSMRb) protein (Recombinant) (STJP000293)

STJP000293

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-31R beta (OSMRb)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse IL-31RB (OSM-R) is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant mouse IL-31 receptor beta, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	18414
Gene Symbol	Osmr
Uniprot ID	OSMR_MOUSE
Immunogen	
Sequence	